
ΔVΔNOS*

GAME READY* GRPRO* 2.1 SYSTEM

User Manual

Systems


- REF 550500-03
- REF 550500-03-DM
- REF 550500-13
- REF 550500-23
- REF 550500-53
- REF 550500-53-DM
- REF 550500-63
- REF 550500-63-DM
- REF 550500-83

Control Units


- REF 550550-03-DM
- REF 550550-03-ED
- REF 550550-03-RN

Description of the GRPro[®] 2.1 System

The following items are included in your GRPro[®] 2.1 System:

A Wrap (comprised of an inner ATX[®] Series Heat Exchanger and an outer Sleeve) must be attached to the system to begin treatment. Each Wrap is sold separately and is not included in the system.

Note: For assistance setting up, using or maintaining the system, or to report unexpected operation or events, please contact Customer Service at 1.888.426.3732 (+1.510.868.2100), from outside of the U.S. please contact your local distributor.

Table of Contents

Description of the GRPro® 2.1 System	EN-2	Storage	EN-8
Introduction To Game Ready®	EN-4	Cleaning.....	EN-8
Register Your Product.....	EN-4	System.....	EN-8
Detailed Instructions for Use	EN-4	Control Unit.....	EN-8
Intended Purpose	EN-4	Connector Hose	EN-8
Intended User	EN-4	Medical Connector Hose.....	EN-8
Patient Target Group.....	EN-4	Wraps	EN-8
Indications For Use.....	EN-4	Carry Bag	EN-8
Clinical Benefits	EN-4	UL Classification.....	EN-9
Specifications	EN-4	Protection Against Electric Shock	EN-9
Pneumatic Compression Cycles.....	EN-4	Pollution Degree Classification.....	EN-9
Temperature Control.....	EN-4	Electromagnetic Interference	EN-9
Contraindications	EN-4	Applied Parts.....	EN-9
Residual Risk and Side Effects.....	EN-5	Electromagnetic Compatibility	EN-9
General Warnings and Cautions.....	EN-5	Troubleshooting	EN-11
Warnings	EN-5	Warranty	EN-13
Cautions.....	EN-5	Manufacturer's Warranty	EN-13
Accessories	EN-5	Warranty Limitations and Exclusions.....	EN-13
Maintenance	EN-6	Limitations of Liability.....	EN-13
User Interface.....	EN-6	Warranty Registration	EN-13
Modes of Operation	EN-6	Registration	EN-13
Buttons.....	EN-6	Extended Warranties	EN-13
Adjusting Temperature	EN-6		
Display.....	EN-6		
Operating the System	EN-7		
Priming the Wrap	EN-7		
Manual Mode	EN-8		
Program Mode	EN-8		

 Medical device	Rx Only	 Caution	 Consult instructions for use	 Dispose of properly
 Serial number	 Catalogue number	 Batch code	 Unique Device Identifier	 Country of manufacture
 Date of manufacture	 Manufacturer	 Importer	 This side up	 Keep dry
 Fragile, handle with care	 33°F 50°C 122°F Temperature limit	 Power plug type	 --- --- --- Direct current 12V/2.5A	 Type BF equipment
 Refer to instructions manual/booklet	 Game Ready® System as to electrical shock, fire and mechanical hazards only in accordance with EN 60601-1:2006 + A1:2013, IEC 60601-1:2012, IEC 60601-1-6:2010(3rd Ed.) + A1:2013, IEC 60601-1-11:2015 and ANSI/AAMI ES60601-1:2005(R)2012.			
IP22 Ingress Protection				

Introduction To Game Ready*

The Game Ready* GRPro* 2.1 System and its accessories, including Wraps, is a medical device that helps athletes and patients recover from post-surgical and acute injuries to reduce edema, swelling, and pain.

Immediately after suffering a musculoskeletal injury, the body initiates a series of physiological responses to defend surrounding tissues and begins to repair the damage. While inflammation is a natural and necessary mechanism in this process, controlling it effectively can actually allow the body to enter the later stages of healing faster. Until now, the RICE (Rest-Ice-Compression-Elevation) principles have been used to passively control symptoms, moderating pain and swelling. Going beyond static cold and compression applications, Game Ready* with ACCEL* Technology proactively aids lymphatic function, encourages cellular oxygen supply, and stimulates tissue repair.

Game Ready* is the ultimate in integrated cold and compression therapy and clinically proven outcomes. The powerful Game Ready* System and complete range of anatomically designed, dual-action wraps allow for optimal therapeutic benefits.

Register Your Product

You must complete your Warranty Registration within 30 days from the purchase date for both the GRPro* 2.1 Control Unit and the Wraps online at www.gameready.com. A Wrap registration card is packaged with each Wrap. Further warranty information can be found in the Warranty Section of this manual.

Don't Just Treat Symptoms, Help Accelerate Healing

We'd enjoy hearing from you. If you have questions, want to learn more, or wish to share your experience with the Game Ready* System, please call the Technical Support Team at **1.888.426.3732**, or email us at GRCustomerService@avanos.com.

1.888.GameReady (1.888.426.3732) | www.gameready.com

Detailed Instructions for Use

Users are required to read this User's Manual in full prior to using this device.

Do not use this device without your physician's specific recommendations for the frequency and duration of your treatments.

The most common recommendation is to use cold therapy at least 4 times a day, for approximately 30 minutes each time, with at least a 30-minute break between treatments.

While the temperature is adjustable, it is reported that the greatest benefit from cold therapy is in the 40-60°F (4.5-15.5°C) temperature range.

Generally reported compression recommendations are for "None" to "Low" pressure settings during the first 24-48 hours after surgery or injury, increasing to "Medium" or "High" only if it is comfortable after the first 48 hours.

We recommend that "High" pressure is never applied if you are confined to bed.

Cautio: United States Federal Law restricts this device to sale by or on the order of a licensed health care practitioner.

Follow the treatment recommendations of your health care practitioner for duration and frequency of use for this device.

Intended Purpose

The Game Ready* GRPro* 2.1 System combines cold and compression therapies with the intention to reduce post-surgical and acute injury-associated pain, edema, and swelling.

Intended User

The GRPro* 2.1 System is intended to be used by or on the order of licensed healthcare professionals in hospitals, outpatient clinics, athletic training settings, or home settings.

Patient Target Group

Patients 18 years and older.

Indications for Use

The GRPro* 2.1 System is intended to treat post-surgical and acute injuries to reduce edema, swelling, and pain for which cold and compression are indicated.

Clinical Benefits


Clinical benefits include reduction of post-surgical- and acute injury-associated pain and edema/swelling.

Specifications

- Size: 16.25" length x 7.75" width x 9.25" height (413 x 197 x 235) mm, not including Carry Bag
- Weight: 7.3 lb. empty, approximately 18 lb. full of ice and water
- AC power: 100-240 V~, 50-60 Hz, 1.6A
- DC input: 12V/2.5 A
- The maximum operating temperature of the equipment is between 33.8-104 °F (1-40 °C).
- The maximum operating altitude of the equipment is 9,843 feet (3,000 meters).
- The expected service life of the Control Unit is 5 years. The expected service life of the parts and accessories shipped with the Control Unit is 2 Years.

Pneumatic Compression Cycles

- High: cyclic 5-75mm Hg
- Med: cyclic 5-50 mm Hg
- Low: cyclic 5-15 mm Hg
- NO Pressure: Wrap shall be vented to atmosphere
- Air Pressure Accuracy: ±10 mm Hg

Temperature Control

The cooling temperature of the circulating ice water will be adjustable between 34 °F (1 °C) and 50 °F (10 °C) as long as the ice water in the reservoir is supplied with sufficient amount of ice.

- Water Temperature Accuracy: ±4 °F (±2.2 °C)

Contraindications

X = Absolute Contraindication – Therapy in these situations should not be used in patients: R = Relative Contraindication – Therapy for these conditions should be used only under the supervision of a licensed healthcare practitioner in patients:	Therapy Mode	
	Cryotherapy	Compression
Who are in the acute stages of inflammatory phlebitis in the affected region.		X
Who have any current clinical signs suggestive of deep vein thrombosis in the affected region.		X
Who have significant arteriosclerosis or other vascular ischemic disease in the affected region.		X
Who have any significant risk factors or current clinical signs of embolism (e.g. pulmonary embolus, cerebral infarction, atrial fibrillation, endocarditis, myocardial infarction, or atheromatous embolic plaque).		X
Who have a condition in which increased venous or lymphatic return is not desired in the affected extremity (e.g., carcinoma).		X
Who have decompensated hypertonia in the affected region.		X
Who have an open wound in the affected region (the wound must be dressed prior to use of Game Ready®).		R
Who have an acute, unstable (untreated) fracture in the affected region.		R
Who are children under 18 years old or patients who have cognitive disabilities or communication barriers, whether temporary (due to medication) or permanent.	R	R
Who have a cardiac insufficiency or congestive heart failure (with associated edema in the extremities or lungs).		R
Who have a localized unstable skin condition (e.g., dermatitis, vein ligation, gangrene, or recent skin graft) in the affected region.	R	R
Who have erysipelas or other active infection in the affected region.		R
Who have significant vascular impairment in the affected region (e.g., from prior frostbite, diabetes, arteriosclerosis or ischemia).	X	
Who have known hematological dyscrasias which affect thrombosis (e.g., paroxysmal cold hemoglobinuria, cryoglobulinemia, sickle-cell disease, serum cold agglutinins).	X	
Who have Raynaud's disease or cold hypersensitivity (cold urticaria).	R	
Who have hypertension or extreme low blood pressure.	R	
Who have diabetes.	R	
Who have compromised local circulation or neurologic impairment (including paralysis or localized compromise due to multiple surgical procedures) in the affected region.	R	
Who have rheumatoid arthritis in the affected region.	R	

Residual Risk and Side Effects

Disclosure of residual risks:

While side effects vary depending on the type of procedure performed, mitigations have been implemented to minimize the following residual risks related to the use of Game Ready® GRPro® 2.1 System:

- Asphyxiation
- Blunt trauma
- Burns
- Choking
- Constricted blood circulation
- Delayed or prolonged treatment / procedure / therapy (no significant medical risk)
- Electric shock
- Hypothermia / frost bite
- Infection
- Inflammatory response
- Interruption of non life-sustaining medical devices
- Interruption of pacemaker or other life-sustaining medical devices
- Laceration
- Malfunction / interruption of other medical devices
- Medical procedure required / secondary surgery
- Nerve palsy / paralysis
- Skin / tissue damage / irritation
- Temporary discomfort

General Warnings and Cautions

Warnings

- Follow the treatment recommendations of your health care practitioner for duration and frequency of use for this device.
- Improper placement or prolonged use of the Game Ready® System could result in tissue damage.
- Use extra caution during the immediate postoperative period, especially when sedated or on any medication that could alter normal pain sensation. Check the skin of the treated region frequently and use mid-to-higher (warmer) temperature range settings or leave more time between treatments, if necessary.
- A layer of clothing between Wrap and skin is recommended for all patients.
- During the course of therapy, patients should monitor the skin of the treated region, the surrounding area and the digits of the extremities of the treated limb (if applicable) for any burning, itching, increased swelling, or pain. If any of these signs are present, or any changes in skin appearance occur (such as blisters, increased redness, discoloration, or other noticeable skin changes), patients are advised to discontinue use and consult a physician.
- Game Ready® Wraps are not sterile; do not place directly against open wounds, sores, rashes, infections, or stitches. The Wrap may be applied over clothing or dressing.
- Game Ready® Wraps are available in multiple configurations but are not intended for all possible physiological uses. For example, the Ankle Wrap is not designed for use on the toes and the Back Wrap is not designed for use in the abdominal region.
- Use extra caution if you are currently on any medication or have an underlying medical condition that predisposes to excessive bleeding or bruising.

Measured Skin Temperature	
Minimum temperature	Skin temperature measured as low as 41 °F (5 °C) when set to coldest therapy

 **WARNING:** To comply with California Proposition 65, the following warning has been included: This product can expose you to chemicals including cadmium, chromium, lead, polybrominated biphenyls or mercury, which are known to the State of California to cause cancer or reproductive toxicity. For more information go to www.p65warnings.ca.gov.

Notification to Users/Patients: Please immediately report any serious incident that may have occurred due to this medical device to the manufacturer or its Authorized Representative and to your National Authority.

Cautions

- To avoid the risk of electrical shock, do not remove any panels from the Control Unit. Opening the case will void the Game Ready® warranty. For all servicing and repair, in the U.S. contact the Technical Support Team at 1.888.426.3732 (+1.510.868.2100) or GRCustomerService@avanos.com; from outside of the U.S. please contact your local distributor.
- To avoid the risk of electrical shock, always turn off the system and disconnect the power line cord from its electrical outlet when not in use or before adding or emptying ice and water.
- Servicing and maintenance shall not be performed while the device is in use.
- Do not use any AC adapter other than that provided by Game Ready. Use of other adapters may result in electrical shock and will void the Game Ready warranty.
- To avoid damage to your product, do not operate the system without any water in the reservoir.
- To avoid electrical shock, product malfunction or damage, never operate the system with damaged power cords or Connector Hoses, or other mechanical damage, or if the unit is otherwise not fully operational.
- To avoid potential damage to your product, do not pour hot water into the reservoir. The system is not designed to operate, and has not been tested, with hot water.
- Do not use anything but ice and water in the reservoir.
- To avoid damage to your product, do not pick up the Control Unit by the lid. Carry the Control Unit using the handle only.
- To avoid potential damage to your product, do not use other manufacturers' wraps with the Game Ready® System.
- To avoid damage to your product, do not operate the Control Unit without a Connector Hose attached.
- To avoid injury, be careful not to trip over the system's power cords and Connector Hose.
- Monitor hose and cord placement when in use. Store hoses and cords when not in use to avoid potential of strangulation.
- Keep the system, including hoses and cords, away from children and pets.
- Use caution when handling small parts to reduce any possibility of swallowing or inhalation
- The GRPro® 2.1 Control Unit is a technical medical device. To avoid damage to your product, handle it with the same care as you would a laptop computer. Do not drop it, kick it or otherwise abuse it unnecessarily. Such abuse will void all Game Ready warranties.
- Do not place the AC Adapter or Battery Pack inside the reservoir for storage or transport.
- No modification of this device is allowed. Any modifications will void the Game Ready warranty.
- There are no special disposal requirements for the Wraps.

This symbol on the Control Unit or its packaging means that this product must not be disposed of with your household waste. To learn where to drop off your electrical and electronic waste, please contact your local city/municipal waste disposal service office or contact Game Ready for assistance.

Accessories

The GRPro® 2.1 Control Unit can be used with any of the following accessories:

- Any Game Ready® Wrap using ATX® Series Heat Exchangers (wraps made by any other manufacturer CANNOT be used with this system)
- Game Ready supplied Mega MDM-030-A120-5396 power supply
- Game Ready supplied Connector Hose
- Game Ready® Carry Bag
- Game Ready® Drain Mode Adapter
- Game Ready® Rechargeable Battery Pack

If long storage of the Battery Pack is required, it is recommended that the battery is periodically recharged to preserve the electronic integrity of the battery. Fully charge before use and every six (6) months.

Maintenance

The reservoir filter should be inspected, cleaned, and/or replaced as necessary.

- 1 Identify the filter within the ice reservoir.
- 2 Using two fingers, grasp and squeeze the two protruding prongs.
- 3 Slide the filter out.
- 4 Rinse debris from the filter and be sure there are no obvious signs of damage. If you have questions or would like to order a new filter, in the U.S. contact the Technical Support Team at 1.888.426.3732 (+1.510.868.2100) or GRCustomerService@avanos.com; from outside of the U.S. please contact your local distributor.
- 5 To replace the filter, or install a new filter, first ensure that the filter is oriented properly with the plastic tab facing up. If the filter is not oriented properly, the protruding tab will prevent the filter from sliding back into place.
- 6 You will feel and hear the filter snap back into place.

Other than maintaining the reservoir filter, no product service or modification should be performed by the user. If your GRPro* 2.1 System requires service, in the U.S. contact the Technical Support Team at 1.888.426.3732 (+1.510.868.2100) or GRCustomerService@avanos.com; from outside of the U.S. please contact your local distributor.

User Interface

Modes of Operation

Manual Mode: The system automatically starts in this mode, and allows the user to adjust treatment time and pressure settings.

Program Mode: This mode allows the user to choose one of six treatment programs that provide therapy for a set time then sleep (no treatment) for a set time, continuously, at a specific pressure setting.

Drain Mode: Allows a user to connect a Hose with a special Hose Adapter (purchased separately) to the unit, enter Drain Mode and have the unit empty the water out of the reservoir through the Hose. Drain Mode can be accessed by pressing the program key and toggling through all six programs. To empty water in the system while in Drain Mode, press the program key until you reach Drain Mode, attach the Hose Adapter to the Connector Hose, place the Hose Adapter over a sink, and press the Play button.

Drain Mode is indicated by the following icons:

Buttons

	Power: Use this button to turn the Control Unit on and off.
	Program: Use this button to select one of the available Programs or to return to Manual Mode. See page 8 in this manual for more information on Programs.
	Play/Pause: Use this button to start or pause a treatment.
	Add Time: Use this button to add time in Manual Mode (does not work in Program Mode). You can add up to 90 minutes. Treatment must be paused in order to add or subtract time.
	Subtract Time: Use this button to reduce time in Manual Mode (does not work in Program Mode). You can subtract up to 90 minutes. Treatment must be paused in order to add or subtract time.
	Pressure Selection: Use this button to select one of four pressure settings: No Pressure, Low Pressure (5-15 mm Hg), Medium Pressure (5-50 mm Hg), and High Pressure (5-75 mm Hg). Pressure selection is not available in Program Mode. Treatment must be paused to change pressure settings.
	Volume: Use this button to select the option of sound or no sound. Push to mute sound. Error notifications will still sound even with Volume off.
	C/F Button: Use this button to select either Celsius or Fahrenheit on the temperature display.
	Backlight: Use this button to turn the backlight on or off.

Adjusting Temperature

To adjust the temperature being applied during treatment, simply turn the temperature knob towards 3 snowflakes for the maximum amount of cold, or towards 1 snowflake for the least amount of cold. Notice that as you adjust the knob, the “Target Temperature” on the display will change. The system will automatically adjust to match the temperature you have selected.

TIP: The maximum amount of cold is dictated by the amount of ice in the reservoir and the setting of the temperature adjustment knob. You may need to stir or replenish the ice to achieve coldest temperatures. Ice will melt faster in the acute phases of injury as the Game Ready* System is removing greater amounts of heat from the treatment site.

Display

Status bar:		
Vol	Vol	Off (Mute) On
Manual Program	Manual Program	Indicates selected mode
Pause	Pause	Indicates unit is paused
Play	Play	Indicates unit is running
Other icons:	Error	Indicates an error. See troubleshooting for error codes
	Sleep	Indicates sleep mode time
	Battery	Indicates remaining battery capacity
Temperature:		Actual Temperature: Approximate temperature of the water leaving the Control Unit
		Target Temperature: Indicates the target temperature setting. The Control Unit will automatically match the Actual Temperature and Target Temperature as closely as possible.
When the system is running in Program Mode:		Displays the amount of run time remaining. Run time count-down bars. Displays the amount of sleep time remaining.
When the system is running in Manual Mode:		Displays the amount of run time remaining.
When selecting a program: 		
Pressure: 		

Operating the System

To operate your GRPro® 2.1 System, you need:

- Control Unit filled with ice and water as indicated by the fill line labels within the reservoir. Optimal performance is achieved by first adding 1.5 liters of water, and then filling the reservoir to the top with ice.
- Game Ready supplied power supply.
- Connector Hose.
- Wrap (includes a Heat Exchanger and Sleeve).

Notes:

- The Wrap is comprised of an inner Heat Exchanger and an outer Sleeve. The combination of Sleeve and Heat Exchanger is referred to throughout this manual as a “Wrap.” To ensure proper performance, be sure to use ATX® Series Heat Exchangers.
- The Control Unit should be placed on a stable surface (such as the floor or a table) during use.
- Note that using the system in an environment with a high ambient temperature may affect its ability to provide adequate cooling, or may limit the ice life.
- If you will be using the system with a Game Ready® optional battery pack (sold separately), please consult the instructions for use that accompany that battery pack.

Priming the Wrap

To prevent 04 errors at the beginning of therapy, it is important to prime the Wrap before a therapy.

- With the GRPro® 2.1 Control Unit off, attach the Connector Hose to the Control Unit and the Wrap
- Lay the Wrap open and flat next to the Control Unit (not on the body)
- Turn the System on and run for 2 minutes with No Pressure
- The wrap is now ready for use and therapy may begin as instructed

Wrap priming instructions may vary depending on Wrap. Refer to the Wrap Use Guide included with each Wrap for specific priming instructions.


Warning:

- When using the Mega MDM-030-A120-5396 power supply, a grounded electrical outlet is not required.
- Position the Control Unit to minimize the risk of tripping over the Control Unit, Connector Hose, or power cord.


Precautions:

- Failure to properly follow the instructions of this manual and those of your medical provider may interfere with or prevent delivery of appropriate therapy.
- To avoid risk of electrical shock unplug the Control Unit from the electrical outlet prior to filling the Control Unit with ice and water.


- 1 Push the door release button to open the reservoir door.

- 2 Add water to fill-line indicated on label within the reservoir. DO NOT OVERFILL. Add ice to top of reservoir.


- 3 Close the reservoir door. Make sure you hear it click.


- 4 Place the Control Unit in the location where you plan to use it. The Control Unit should only be used in the upright orientation as shown. The Control Unit will leak if placed on its side.

- 5 Connect the AC Adapter to the receptacle on the end panel of the Control Unit, then plug the AC Adapter into an electrical outlet. The power indicator light (orange color) on the Control Unit should illuminate. Press the power button. The screen should light up and the Control Unit should beep twice. The power indicator light should turn from orange to green. If you do not see these indicators, refer to “Control Unit will not turn on” on page 12 of the Troubleshooting Guide in this manual.

Note: Make sure not to position the device so that it is difficult to operate the plug.

- 6 Connect the larger end of the Connector Hose (with the red button) to the Control Unit. Make sure you hear it click. To disconnect, simply press the red button and remove the connector from the Control Unit. Note: The Control Unit should be turned on prior to attaching a Wrap.


- 7 Apply the selected Wrap (consult the Use Guide accompanying each Wrap for application instructions). To ensure proper performance of the system, it is important to use ATX® Series Heat Exchangers in your Wrap.

Caution:

- The Wrap is not sterile. Do not place directly against open wounds, sores, rashes, infections, or stitches. The Wrap may be applied over clothing or dressing.
- To ensure best fit, be sure the Wrap is completely deflated prior to each application.

- 8 Connect the smaller end of the Connector Hose (with a blue or gray button) to the Wrap. Make sure you hear it click. To disconnect, simply press the blue or gray button and remove the connector from the Wrap.

If you have any problems with the set-up of your GRPro® 2.1 System, in the U.S. contact the Technical Support Team at 1.888.426.3732 (+1.510.868.2100) or GRCustomerService@avanos.com; from outside of the U.S. please contact your local distributor.

Warnings:

- Follow the treatment recommendations of your health care practitioner for the use of this device.
- Improper placement or prolonged use of the GRPro® 2.1 could result in tissue damage.
- READ COMPLETE INDICATIONS, CONTRAINDICATIONS, CAUTIONS AND WARNINGS BEFORE USING THIS PRODUCT!

9 Manual ▶ **6** Select the operating mode. You can choose either Manual Mode or Program Mode. Manual Mode allows you to customize the time and pressure settings. Program Mode allows you to select one of six automated programs that provide treatment for a set time interval then “sleep” (providing no treatment) for a set time, continuously, at a specific pressure setting (refer to list of available programs below). The unit automatically starts in Manual Mode.

Manual Mode

 	Set the time in five minute increments by pushing the +/- buttons.
	Set the pressure by pushing the pressure button. You can select from 4 pressure settings: No pressure, Low Pressure (5-15mmHg), (Medium Pressure (5-50mmHg), High Pressure (5-75mmHg).

Program Mode

	Push the program button to enter Program Mode. In Program Mode, the unit will operate continuously according to the selected program. You will need to replenish ice and water as previously indicated in step 2.
	You can select from the following programs: (Push the Program Button to scroll through the available programs.) Program 30: 30 minutes on, 30 minutes sleep. No pressure. Program 2: 30 minutes on, 30 minutes sleep. Low pressure. Program 3: 30 minutes on, 30 minutes sleep. Medium pressure. Program 4: 30 minutes on, 60 minutes sleep. No pressure. Program 5: 30 minutes on, 60 minutes sleep. Low pressure. Program 6: 30 minutes on, 60 minutes sleep. Medium pressure. Program d: Drain Mode. Please refer to page 6 in this manual for details.

10  Press play/pause button to start your Game Ready* treatment with ACCEL* Technology. Press play/pause at any time to stop your treatment.

11 Set to your target temperature (displayed in the Target Temperature window) by turning the knob. For coldest therapy, turn the knob clockwise towards three snowflakes. For warmer therapy, turn the knob counter-clockwise towards one snowflake.

12 Replenish ice and water levels as necessary, and indicated in Step 2, to maintain your target temperature.

13 After use:

- Unplug the AC Adapter and the Connector Hose from the Control Unit
- Push the door release button to open the door
- Carefully pour out the ice and water
- Wipe off any excess water

Storage

When you are done using the system for the day:

- Unplug the AC Adapter and the Connector Hose from the Control Unit.
- Push the door release button to open the door.
- Carefully pour out the ice and water.
- Wipe off any excess water.
- Store the Control Unit with the lid fully open to allow the interior to dry and preserve the reservoir seal. Store the Control Unit in the Carry Bag or in another safe location. Remember that your GRPro* 2.1 Control Unit is a valuable piece of equipment and should be treated with great care, like a laptop computer.

System Storage Temperature: 33 °F - 122 °F (1 °C - 50 °C)

Relative Humidity: 15% - 90% non-condensing

Caution: Do not keep in extreme cold or hot temperatures (below 33 °F or above 122 °F or below 1 °C or above 50 °C). Do not leave in a hot or freezing car. Do not leave the Control Unit in direct sunlight. The UV light may damage or discolor the Control Unit.

Cleaning

System

If desired, the system (Control Unit, Hoses, and Wraps) can be flushed with a mixture of isopropyl alcohol and water, and run on no compression (“off”). This should be followed by running the Control Unit with clean water only on the no compression (“off”) setting.

Control Unit

The exterior of the Control Unit and the visible interior surfaces of the reservoir can be cleaned with a soft cloth and one of the following cleaning agents:

- Mild detergent
- 70% Isopropyl alcohol
- Quaternary ammonium based cleaning agents (such as Cavicide[®])

Procedure:

- Follow the manufacturer’s instructions and precautions for the cleaning agent you select.
- Apply the selected cleaner to a soft cloth and wipe down all surfaces of the Control Unit.
- Allow the Control Unit to dry thoroughly before storing it in the bag.
- The Control Unit should be cleaned as needed.

DO NOT USE:

- Phenolic-based disinfectants (such as Amphyll[®] – typically only found in a clinical use setting). Doing so will damage the plastics and will void your warranty.
- Any solvent-based cleaners on the Control Unit. Doing so will damage the plastics and will void your warranty.
- Abrasive materials to clean the Control Unit. Doing so will damage the plastics and will void your warranty.

Caution: The Control Unit is not a waterproof device. Do not apply a direct stream of any liquid onto the Control Unit, submerge the Control Unit, or allow any liquid to pool on the surface of the front panel of the Control Unit.

Connector Hose

The surface of the Connector Hose can be cleaned using a soft cloth and one of the following:

- Mild detergent
- Steri-Fab[®]
- 70% Isopropyl alcohol

DO NOT USE:

- Quaternary ammonium based cleaning agents (such as Cavicide[®]). Doing so will damage the plastics and will void your warranty.
- Phenolic-based disinfectants (such as Amphyll[®] – typically only found in a clinical use setting). Doing so will damage the plastics and will void your warranty.
- Any solvent based cleaners. Doing so will damage the plastics and will void your warranty.
- Abrasive materials. Doing so will damage the plastics and will void your warranty.
- Any petroleum based lubricants. Doing so will damage the o-rings and will void your warranty. If lubrication is required, the use of silicon spray is recommended.

Medical Connector Hose

(550500-13 and 550500-63 only)

The surface of your Game Ready* Medical Connector Hose may be wiped down as needed using soap and water.

To disinfect the surface, certain commercially available agents may be used, always following the instructions provided with the agent. We recommend Quaternary Ammonia based cleaning agents, such as Cavicide[®]. When possible, we recommend the use of Caviwipes[®]. A liquid form may also be applied using a soft cloth. In all cases, follow the cleaning agent manufacturer’s instructions to determine the necessary exposure time for proper disinfection.

Wraps

Gently remove Heat Exchanger from the Sleeve. Hand wash the Sleeve in cold water, using a mild detergent or antibacterial soap. Do not use fabric softener. Hang to dry.

If needed, the external surface of the Heat Exchanger may be cleaned by wiping down with commercial non-bleach cleaning wipes or hand washed using a very small amount of mild detergent or antibacterial soap. **DO NOT MACHINE WASH.**

Refer to the Wrap Use Guide accompanying individual Wraps for more information.

Carry Bag

The Carry Bag should be cleaned using a soft cloth or brush and a mild detergent. Febreze[®] or the equivalent can be used on the Bag if desired. If the Carry Bag has a biological material on the surface, Steri-Fab[®] may be used to disinfect those surfaces.

Be sure to test any product on a small portion of the Bag to ensure that it will not cause damage.

Note: To operate the GRPro* 2.1 System, you do not need to remove it from the Carry Bag. Simply unzip the Bag’s main compartment and end panel. Fill the reservoir with ice and water. Attach the Connector Hose and the AC Adapter to the end panel of the Control Unit and plug the AC Adapter into an electrical outlet.

UL Classification

Protection against electric shock:

The GRPro* 2.1 System is considered to be Class II when connected to the Mega MDM-030-A120-5396 power supply.

The MDM-030-A120-5396 power supply is intended for home use.

Protection against harmful ingress of water:

This product provides ordinary protection against ingress of water. When used with the MDM-030-A120-5396 power supply the device is classified to an IP22 rating, which is a protection against dripping water.

Pollution degree classification:

This product is classified as Pollution degree 2.

Degree of safety in the presence of flammable anesthetics or oxygen:

Not suitable for use in an oxygen enriched environment or in the presence of flammable anesthetics.

Electromagnetic interference:

This equipment has been tested and found to comply with the limits for medical devices in IEC 60601-1-2 Ed. 4.1 b:2020. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other devices

in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to other devices, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving device.
- Increase the separation between the equipment. Connect the equipment into an outlet on a circuit different from that to which the other device(s) are connected.
- Consult the manufacturer or field service technician for help.

The device is intended for use in orthopedic centers, therapy clinics, athletic training facilities, hospitals, nursing facilities, medical centers, as well as in patient homes. The device should not be used in an environment where the intensity of electromagnetic disturbance is high.

In the event the device experiences a loss of performance or degradation due to electromagnetic disturbance, the device is expected to continue to operate safely.

The device should not be used less than 30 cm (12 inches) away from any portable and mobile RF communications equipment.

Applied Parts:

The following are classified as Type BF Applied Part:

- Patient Wrap
- Hose and Connectors
- Control panel (Membrane Buttons + Knob)^{††}

^{††} Parts may reach up to 119 °F (48 °C). Sustained contact duration not to exceed 1 minute.

Electromagnetic Compatibility

Electromagnetic Compatibility Information According to IEC/EN 60601-1-2

Table 1 for Emissions

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS		
The GRPro* 2.1 is intended for use in the electromagnetic environment specified below. The customer or the user of the GRPro* 2.1 should assure that it is used in such an environment.		
EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
RF emissions CISPR 11	Group 1	The GRPro* 2.1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The GRPro* 2.1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Table 2 for Electromagnetic Immunity

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY			
The GRPro* 2.1 is intended for use in the electromagnetic environment specified below. The customer or the user of the GRPro* 2.1 should assure that it is used in such an environment.			
IMMUNITY TEST	PROFESSIONAL HEALTHCARE FACILITY ENVIRONMENT LEVEL	HOME HEALTHCARE ENVIRONMENT LEVEL	GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact	± 8 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
	± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Electrical fast transient/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	N/A
Surge IEC 61000-4-5 ±0.5kV, ± 1kV, ± 2kV line(s) to ground	±0.5kV, ± 1kV, line(s) to line(s)	±0.5kV, ± 1kV, line(s) to line(s)	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% <i>U_i</i> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% <i>U_i</i> ; 1 cycle and 70% <i>U_i</i> ; 25/30 cycles Single phase: at 0° 0% <i>U_i</i> ; 250/300 cycle	0% <i>U_i</i> ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% <i>U_i</i> ; 1 cycle and 70% <i>U_i</i> ; 25/30 cycles Single phase: at 0° 0% <i>U_i</i> ; 250/300 cycle	If the user of the GRPro* 2.1 requires continued operation during power mains interruptions, it is recommended that the GRPro* 2.1 be powered from an uninterruptible power supply or a battery.
Rated Power Frequency Magnetic Fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz	N/A
NOTE UT is the AC mains voltage prior to application of the test level.			

Table 3 for RF Electromagnetic Immunity

GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY			
The GRPro* 2.1 is intended for use in the electromagnetic environment specified below. The customer or the user of the GRPro* 2.1 should assure that it is used in such an environment.			
IMMUNITY TEST	PROFESSIONAL HEALTHCARE FACILITY ENVIRONMENT LEVEL	HOME HEALTHCARE ENVIRONMENT LEVEL	GUIDANCE
Conducted Disturbances Induced by RF Fields IEC 61000-4-6	3V 0.15 MHz - 80 MHz 6V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz See Table 4	3V 0.15 MHz - 80 MHz 6V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz See Table 4	Portable and mobile RF communications equipment should be used no closer to any part of the GRPro* 2.1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz
Radiated RF EM Fields IEC 61000-4-3	3 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	10 V/m 80 MHz - 2.7 GHz 80% AM at 1 kHz	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).
Proximity Magnetic Fields IEC 61000-4-39	134.2 KHz, Pulse Modulation 2.1 kHz, 65 A/m 13.56 MHz, Pulse Modulation 50 kHz, 7.5 A/m	30 kHz, CW Modulation, 8 A/m 134.2 KHz, Pulse Modulation 2.1 kHz, 65 A/m 13.56 MHz, Pulse Modulation 50 kHz, 7.5 A/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the GRPro* 2.1 is used exceeds the applicable RF compliance level above, the GRPro* 2.1 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the GRPro* 2.1. ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 4 for Immunity to RF Wireless Communications Equipment

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	IMMUNITY TEST LEVEL (V/m)
385	380 to 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	27
450	430 to 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	28
710	704 to 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	9
745				
780				
810	800 to 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	28
870				
930				
1,720	1,700 to 1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	28
1,845				
1,970				
2,450	2,400 to 2,570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	28
5,240	5,100 to 5,800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	9
5,500				
5,785				
If necessary to achieve the IMMUNITY TEST LEVEL, the distance between transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.				
^{a)} For some services, only the uplink frequencies are included. ^{b)} The carrier shall be modulated using a 50% duty cycle square wave signal. ^{c)} As an alternative to FM modulation, the carrier may be pulse modulated using a 50% duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case.				

Table 5 for RF Immunity

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE GRPRO* 2.1			
The GRPro* 2.1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the GRPro* 2.1 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the GRPro* 2.1 as recommended below, according to the maximum output power of the communications equipment.			
RATED MAXIMUM OUTPUT POWER OF TRANSMITTER W	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER		
	150 KHZ TO 80 MHZ D = 1.2 √P	80 MHZ TO 800 MHZ D = 1.2 √P	800 MHZ TO 2.5 GHZ D = 2.3 √P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

WARNINGS:

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the User's Manual.

Portable and mobile RF communications equipment can affect Medical Electrical Equipment.

The GRPro* 2.1 System has to be powered with the AC adapter Mega MDM-030-A120-5396 power supply in order to be compliant with IEC/EN 60601-1-2 section 6.1 and 6.2.

The use of accessories, transducers and cables other than those specified and sold by the manufacturer of the GRPro* 2.1 System as replacement parts for internal components may result in increased emissions or decreased immunity of the medical electrical system.

The GRPro* 2.1 System should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the GRPro* 2.1 System should be observed to verify normal operation in the configuration in which it will be used.

Troubleshooting

ERROR	WHAT DOES IT MEAN?	WHAT CAN I DO?
 01	Air Pressure Sensor: The Control Unit has detected a problem calibrating the air pressure circuit on start-up.	• This is most likely to occur if you have restarted the system with an inflated Wrap attached. • Detach the Wrap, pressing it flat to expel the air accumulated inside and try again.
 02	Self-Test Error – Air Pump: The Control Unit has detected an electronic problem in the air circuit on start-up.	• Disconnect the Wrap from the Control Unit. • Turn the Control Unit off and on again without a Wrap connected. • Reconnect the Wrap and resume treatment. • If the problem persists, contact The Technical Support Team.
 04	Dry Pump: The Control Unit has detected a dry pump. In order to prevent possible damage to the fluid pump, the unit will stop therapy.	• Be sure you are using ATX* Series Heat Exchangers. • Note that if you are using new Wraps, the water in the reservoir may have been depleted and causing this error. Make sure there is adequate water in the reservoir based on the fill line indicator label. Refer to instructions for priming a Wrap below. • Verify that the reservoir filter is not clogged (refer to filter maintenance instructions on page 6). • Make sure there are no kinks in the Wrap or Connector Hose. • Turn the Control Unit off and on again. • Disconnect and reconnect the hose from the Control Unit and the Wrap, verifying that an audible “click” is heard at both connection points. • Prime the Wrap using the following steps: • Select “Off” (No Pressure). • Attach the hose to the unit and the Wrap. • Lay the Wrap open and flat next to or lower than the Control Unit (not on the body). • Run the system for 2 minutes. • Prime the Control Unit using the following steps: • Disconnect the hose from the Control Unit. • Now, look at the Wrap connection location on the Control Unit. On the top valve, push the white prong in so that it is flush with the metal connector. • Make sure you are not fully covering the opening in the prong. • Push start, and water should squirt out of the valve. • Restart the system.
 06	Over Pressure: Indicates that the Control Unit has exceeded the target air pressure.	• Turn the Control Unit off and back on. • Make sure the Wrap is applied securely. • Do not make sudden movements during treatments. Rapid shifting in position may produce a quick change in pressure in the Wrap and cause this error.

 07	Under Pressure: The Control Unit can't reach its target maximum compression. This often indicates that there is a leak in the pneumatic compression circuit, either in the Connector Hose, Wrap or Control Unit. Or it may occur because the hook and loop fastener on your Wrap has worn out.	• Make sure the Wrap is applied securely. • Try using a different Wrap and Hose to isolate which component may be producing the error. For example, an error which occurs with a Shoulder Wrap but not a Knee Wrap may indicate that the Shoulder Wrap is causing the error: not the Control Unit. • If using a Dual Hose, make sure you have two Wraps attached.
 08	Deflation Error: The Control Unit has detected that the Wrap has not properly deflated.	• Turn the Control Unit off. • If there is air left in the Wrap, disconnect the Wrap from the Connector Hose and manually deflate the Wrap by applying pressure to it. • Reconnect the Hose to the Wrap and reapply the Wrap to the body. • Turn the Control Unit on and press play/pause. • Make sure the Wrap is applied securely against the body. • If possible, try using a different Wrap and Hose to isolate which component may be producing the error. For example, an error which occurs with a Shoulder Wrap but not a Knee Wrap may indicate that the Shoulder Wrap is causing the error: not the Knee Wrap or Control Unit.
 09	Pump Performance Error: The Control Unit has determined that the fluid pump may be working too hard. This could be caused by ice or debris in the fluid circuit. In order to prevent possible damage to the fluid pump, the unit will stop therapy.	• Turn the Control Unit off and back on again. • Reapply the Wrap, making sure to follow all application instructions accompanying the Wrap. • Disconnect and reconnect the hose from the Control Unit and the Wrap, verifying that an audible “click” is heard at both connection points. • If that does not solve the problem turn the Control Unit off for 20 minutes (to let the pump cool down) before turning it on again to try again.
 10	Low Flow: Control Unit has detected something blocking the water flow.	• Be sure you are using ATX® Series Heat Exchangers. • Check all hose connections. • Disconnect and reconnect Wrap from Connector Hose. • Make sure there is water in the reservoir. • Verify that the reservoir filter is not clogged. • Make sure there are no kinks in the Wrap or Connector Hose. • Reapply the Wrap snugly, making sure to follow all application instructions accompanying the Wrap. • Turn the Control Unit off and on again. • Disconnect and reconnect the hose from the Control Unit and the Wrap.
 12	Self-Test Error – Fluid Circuit: The Control Unit has detected an electronic problem in the fluid circuit on start-up.	• Turn off the Control Unit. • Turn the Control Unit back on and resume treatment. • If the problem persists, contact the Technical Support Team
 13	Temperature Calibration Error: The Control Unit has detected a failure in the Temperature Control circuit or is being operated outside of the advised temperature range (33.8-104 degrees Fahrenheit).	• Be sure you are operating the system within the advised temperature range of 33.8-104 degrees Fahrenheit and storing the system within the advised temperature range of 33-120 degrees Fahrenheit. • Turn the Control Unit off. • Fill the Control Unit with ice and water per the label instructions inside the reservoir. • Turn the unit back on again and push Play/Pause. • Repeat this process up to three times. • If the problem persists, contact the Technical Support Team.
	Control Unit will not turn on:	• Press the power button. If no orange or green light is illuminated, make sure the AC Adapter is plugged into a working electrical outlet and all connections are secure. (Test outlet with another electrical device.) • Press the power button again. If the Control Unit still will not turn on, contact the Technical Support Team.
	Control Unit won't reach Target Temperature, or temperature is unstable:	• The fill-line labels within the reservoir indicate an ice/water ratio that will assist the Control Unit in achieving the temperature you've specified using the Temperature Knob. If you have filled the Control Unit according to the fill-line labels and are still unable to achieve your desired target temperature, try the following steps: • If you are still unable to achieve WARMER temperatures, be sure the Temperature Knob is turned to full warm, use less ice, and reduce the amount of water if necessary. • If you are still unable to achieve COLDER temperatures, be sure the you are using ATX® Series Heat Exchangers. Confirm that the Temperature Knob is set to full cold and that the reservoir is full of ice, replenishing frequently and stirring the reservoir to break up large formations of ice if necessary. Make sure there are no kinks in the Wrap or Connector Hose. Reapply the Wrap, making sure to follow all accompanying application instructions. Finally, if you still cannot achieve the coldest temperature desired, adding more water than what is indicated by the reservoir fill-line labels may also allow the Control Unit to reach colder temperatures. To achieve this effect, water may be added to the top of the ice level. CAUTION: By overfilling the reservoir as in the last step, the temperature control feature of the system will not work and the Control Unit will be administering therapy at “full cold.” Consult the Warnings on page 5 to make sure you take proper steps to minimize the risk of injury.

Need more help with a problem? In the U.S. contact the Technical Support Team at 1.888.426.3732 (+1.510.868.2100) or GRCustomerService@avanos.com; from outside of the U.S. please contact your local distributor.

Warranty

Manufacturer’s Warranty

Avanos Medical, Inc. warrants that the GRPro* 2.1 Control Unit, if properly used, will be free from defects in material and workmanship for a period of one (1) year after the date the GRPro* 2.1 Control Unit was purchased, unless otherwise stipulated by local law or regulation. If the GRPro* 2.1 Control Unit, which is the subject of this Limited Warranty, malfunctions during the warranty period for reasons covered by this Limited Warranty, Avanos Medical, Inc., at its options, will:

- REPAIR the GRPro* 2.1 Control Unit OR
- REPLACE the GRPro* 2.1 Control Unit with another GRPro* 2.1 Control Unit.

THIS LIMITED WARRANTY AND ANY IMPLIED WARRANTIES THAT MAY EXIST UNDER STATE LAW APPLY ONLY TO THE ORIGINAL PURCHASER OF THE GRPRO* 2.1 CONTROL UNIT AND ARE NON-TRANSFERABLE.

Extent of Limited Warranty

This limited warranty does not cover damages due to external causes, including, without limitation, accident, usage not in accordance with product instructions, misuse, neglect, alteration or repair.

How to Obtain Warranty Service

To obtain warranty service, in the U.S. call the Technical Support Team at 1.888.426.3732 (+1.510.868.2100); from outside of the U.S. please contact your local distributor. If you are not sure who the local distributor is, contact the Technical Support Team at +1.510.868.2100 and we will refer you to one. You must have returned the Warranty Registration card to Avanos Medical, Inc. within thirty (30) days from the date of purchase to qualify for warranty service. If you qualify for warranty service from Avanos Medical, Inc., you will be issued a Returned Material Authorization (RMA) number. When you return the GRPro* 2.1 Control Unit to Avanos Medical, Inc., you must write the RMA number on the outside of the package. Avanos Medical, Inc. will not accept returned GRPro* 2.1 Control Units without an RMA number on the package. If you return the GRPro* 2.1 Control Unit to Avanos Medical, Inc., you must assume the risk of damage or loss during shipping. You must use the original packaging or the equivalent. Avanos Medical, Inc. may require you to verify in writing that you are the original purchaser of the GRPro* 2.1 Control Unit. Avanos Medical, Inc. may elect to replace or repair the GRPro* 2.1 Control Unit with either a new or reconditioned product. The returned product shall become Avanos Medical, Inc.’property upon receipt. The replacement GRPro* 2.1 Control Unit is warranted under this written warranty and is subject to the same limitations and exclusions for the remainder of the original warranty period. THIS WARRANTY IS NOT TRANSFERABLE.

Warranty Limitations and Exclusions

THESE WARRANTIES REPLACE ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. AVANOS MEDICAL, INC. MAKES NO EXPRESS WARRANTIES BEYOND THOSE STATED HERE. AVANOS MEDICAL, INC. DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING, WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OF IMPLIED WARRANTIES SO THIS LIMITATION MAY NOT APPLY TO YOU. ALL EXPRESS AND IMPLIED WARRANTIES ARE LIMITED IN DURATION TO THE LIMITED WARRANTY PERIOD. NO WARRANTIES APPLY AFTER THAT PERIOD. SOME JURISDICTIONS DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THIS LIMITATION MAY NOT APPLY TO YOU.

Limitations of Liability

AVANOS MEDICAL, INC.’ RESPONSIBILITY UNDER THIS, OR ANY OTHER WARRANTY, IMPLIED OR EXPRESS, IS LIMITED TO REPAIR OR REPLACEMENT, AS SET FORTH ABOVE. THESE REMEDIES ARE THE SOLE AND EXCLUSIVE REMEDIES FOR ANY BREACH OF WARRANTY. AVANOS MEDICAL, INC. IS NOT RESPONSIBLE FOR DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY BREACH OF WARRANTY OR UNDER ANY OTHER LEGAL THEORY INCLUDING, BUT NOT LIMITED TO, LOST PROFITS, DOWNTIME, GOODWILL, AND DAMAGE TO OR REPLACEMENT OF EQUIPMENT AND PROPERTY. SOME JURISDICTIONS DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATIONS OR EXCLUSIONS MAY NOT APPLY TO YOU. THIS LIMITED WARRANTY GIVES YOU SPECIFIC RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS THAT VARY FROM JURISDICTION TO JURISDICTION.

Warranty Registration

The GRPro* 2.1 Control Unit comes with a 1-year warranty from date of purchase, unless otherwise stipulated by local law or regulation. The Heat Exchanger, Connector Hose, AC Adapter, and Power Cord come with a 1-year warranty. In the case of a manufacturer’s defect, the Wrap sleeve and sleeve accessories may be returned within 7 days of purchase.

Registration

Please complete the Warranty Registration within 30 days for both the GRPro* 2.1 Control Unit and the Wraps online at www.gameready.com. You will need the following information to complete your Warranty Registration: The Control Unit’s model number (REF) and its serial number (SN). These numbers are located on the label on the bottom of the Control Unit. Simply go to www.gameready.com, visit the Product Registration page, fill out the form and submit your information.

Extended Warranties

Extended Warranties are available for the GRPro* 2.1 System. For details and information, in the U.S. call the Technical Support Team at 1.888.426.3732 (+1.510.868.2100); from outside of the U.S. please contact your local distributor.

The CE mark is valid only if it is also printed on the device label.

Distributed in the USA by Avanos Medical Sales, LLC, 5405 Windward Parkway, Alpharetta, GA 30004 USA.
In USA, 1-888-GAMEREADY (1-888-426-3732). www.gameready.com

 Avanos Medical, Inc.,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.

 Avanos Medical Belgium BV,
Leonardo da Vincilaan 1, 1930
Zaventem, Belgium.

Sponsored in Australia by Avanos Medical
Australia Pty Ltd, 475 Victoria Avenue,
Chatswood, NSW 2067, Australia.

eifu.avanos.com
1-844-428-2667

[†] Subject to third party ownership of trademark rights.

*Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates.
© 2018-2024 AVNS. All rights reserved. 2024-10-08 15-M1-2024-02

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The interface shows a large central circular display with a pressure gauge icon and a temperature display. Surrounding it are various control buttons and indicators. Labels point to specific features:

- Pressure (Πίεση):** A gauge with settings: off, Low, Med, High.
- Volume:** A speaker icon with a volume slider.
- Manual:** A button with a '6' icon.
- Program:** A button with a '6' icon.
- Pause:** A button with a pause icon.
- Play:** A button with a play icon.
- Sleep Time:** A button with a smiley face icon.
- Run Time:** A button with a smiley face icon.
- Thermometer:** A vertical thermometer showing temperature in °F/°C.
- Target Temperature:** A display showing 88°F.
- Pressure (Πίεση):** A gauge with settings: off, Low, Med, High.
- Volume:** A speaker icon with a volume slider.
- Manual:** A button with a '6' icon.
- Program:** A button with a '6' icon.
- Pause:** A button with a pause icon.
- Play:** A button with a play icon.
- Sleep Time:** A button with a smiley face icon.
- Run Time:** A button with a smiley face icon.
- Thermometer:** A vertical thermometer showing temperature in °F/°C.
- Target Temperature:** A display showing 88°F.

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