

EverTrac ET-800 Traction unit

The ET-800 is a digital traction unit for continuous, intermittent or harmonic lumbar and cervical traction therapy. It is a no-nonsense machine with all must-have functionality and rock solid build quality. Together with the great user friendliness, high safety standards and nice design the ET-800 will meet your highest professional requirements.



Standard features

- Three therapy modes: continuous, intermittent and harmonic intermittent traction.
- Automatic tension compensation.
- Maximum traction force 91 kg for lumbar traction 20 kg for cervical traction.
- Step less adjustable traction speed.
- Timer from 1 -99 minutes, with remaining time displayed.
- Traction force displayed in kg or lb.
- Digital display for preset force and actual force.
- Digital displays for treatment phase (progressive - traction - regressive - pause).
- Buzzer with adjustable sound level for end of treatment or protection alarms.

Specifications

Dimensions (w x l x h)	: ± 31 x 28 x 26 cm
Weight	: ± 11.4 kg
Mains voltage	: 220 - 240 VAC 50 Hz, 0.6 A
Power consumption	: Max. 150 W
Insulation class	: I type BF (IEC 60601-1)
EMC tested	: IEC 60601-1-2
Classification:	: Class IIa MDD 93/42/EEC Annex IX
CE marking	: CE 2460

Ordering data

160.985 EverTrac ET-800 Traction unit
Standard accessories:
Mains lead (straight, Europe)
Emergency stop switch
Extension cord with hooks
Four knobs for fixation to the Extension frame
Treatment Guide
User manual

Optional accessories

100.457 Neck harness according to Glisson
100.458 Transverse suspension bar
100.459 Set of washable pads
100.476 Hip belt with Velcro fastening
100.477 Chest belt with Velcro fastening
110.487 Flexion stool
160.996 G2 Traction Extension frame for ET-800



100.457 + 100.458 + 100.459



100.476 + 100.477



110.487



160.996 + 160.985