

gymna[®]

SPECIAL EDITION

gymna.GO **Carbon Black**

Basically perfect. Now Boldly Black

TECHNISCHE DATASHEET

Technische specificaties

Afmetingen (l x b)	205 x 67 cm
Gewicht	+/- 78 kg (D1) / +/- 79 kg (T3x) / +/- 81 kg (T5x)
Hoogteverstelling	50 - 100 cm
Hefvermogen	Max. 200 kg
Motorkracht	6 kN
Motorgebruik	Max. 2 min. on / 18 min. off
Netspanning	100 - 230 VAC, 50 - 60 Hz
Stroomverbruik	Max. 2,0 A
Isolatieclassificatie	Class II
Toegepaste onderdelen	Type BF
Vochtbestendig	IP 21
Frame	Stalen frame, gepoedercoat RAL 9004 (zwart)
Kussenafwerking	Lineair, enkelvoudig gewatteerd 40 kg/m ³
Bekleding	Premium PU-inkt zwart Ftalatenvrij, biocompatibel, slijtvast en vlamvertragend.
Classificatie toestel	Class I
Conformiteit 	IEC/EN 60601-1 IEC/EN 60601-1-12 EG-richtlijnen 2017/745 (MDR) & 2012/19/EG 2011/65/EU
Garantie	2 jaar

gymna.GO Carbon Black


functionaliteiten

gymna.GO Carbon Black Series	D1	T3x	T5x
2-delige tafel	●		
3-delige tafel		●	●
Armsteunen			●
i-Control elektrische hoogteverstelling	●	●	●
Veiligheidsschakelaar	●	●	●
Ergomax gezichtsopening	●	●	●
Wielset met vergrendelingsmechanisme	●	●	●
Hoogwaardige PU-bekleding	●	●	●




Behandeltafel afmetingen

gymna.GO D1




Behandeltafel afmetingen

gymna.GO T3x



Behandeltafel afmetingen

gymna.GO T5x



MECO010 - AFP001 - 825.071 - V.100 - 2025 GYMNA.GO CBE TECHNICAL DATASHEET - NL © GYMNAUNIPHY NV