

RCM SERIES

RCM Pro MED G6515BH

BH

P.V.P.

Treadmills

The RCM Pro MED treadmill has an optimised AC motor that reaches speeds from 0.3 to 18km/h, it has been designed for semi-professional use. The side bars make it perfect for rehabilitation purposes.



Monitor

Measures:

Time, inclination, heart rate, speed, distance, calories

- Max. user weight: 180Kg
- Weight: 125Kg
- Dimensions: 185cm x 82cm x 154cm
- Power: 3,5 CV-AC
- Speed: 0.3-18Km/h
- Running surface: 155x55cm



3.5 HP AC ENGINE

3.5 hp AC motor that reaches a maximum speed of 0.3 to 18 km/h, allowing efficient workouts adapted to different levels of intensity and physical condition.



COMFORT BAND

Increases walking comfort by reducing the impact on the sole of the foot.



SIDEBARS

The treadmill includes auxiliary side bars to assist the user during the rehabilitation process.



ELECTRICAL INCLINATION

Adjust the incline level up to 15% for a more intense workout.



Specs	RCM Pro MED G6515BH
Use frequency	Semi-professional
Maximum user weight	180Kg
Power (Peak/Continuous)	3.5 HP-AC
Speed	0.3-18Km/h
Max. electric inclination	15%
ECO mode	Yes
Instant Speed Keys	Yes
Inclination Instant Keys	Yes
Running surface (L x W)	155x55cm
Damping	Yes
Contact pulse measurement	Yes
Fan	No
Soft Drop System (SDS)	No
Transport wheels	Yes
Length	185cm
Width	82cm
Width	154cm
Weight	125Kg
Foldable	No
Programs	
Preset Programs (Prg)	36
Random Program (RP)	No
Customizable profiles (uprg)	3
Fitness Test (FT)	No
Heart rate control program (HRC)	3
Recovery test (RT)	No
Body fat test (BF)	Yes
Monitor	
Screen	LED Dot matrix
Telemetric pulse measurement	Yes, optional band
Bluetooth Pulse	Yes, optional band
iConcept	Yes

Notes